

Rat Mitochondrial DNA Copy Number

Assay Kit

(44 reactions)

Catalog Number: MCN 2

Store at -20°C.

FOR RESEARCH USE ONLY



Introduction: This DNA analysis kit is for the determination of rat mitochondrial DNA copy number, *in vivo* and *in vitro*, by the comparison of mitochondrial (mt) and nuclear (n) DNA measured by real-time PCR.

Kit Contents:

- 96 well PCR plate
- rtPCR reaction mix.
- Validated primers to quantify mitochondrial DNA (mtDNA).
- Validated primers to quantify nuclear DNA (nDNA).
- Positive control [1.825 ng/ μ l] (isolated total DNA from rat NRK-52E cells).

Not Included in Kit:

- DNA isolation Kit
- Nuclease-free water
- PCR Tubes and Caps

Thermal cycler program:

- Preprogram PCR machine for this profile:
 - a. 95°C, 10 min
(40 Cycles)
 - b. 95°C, 15 sec
 - c. 60°C, 60 sec

Real time PCR procedure: The following procedure is for each 20 μ L reaction. Increase all amounts proportionally according to the total number of tubes.

- Per PCR tube (20 μ L Rx), mix the following:
 - a. 1 μ L forward primer
 - b. 1 μ L reverse primer
 - c. 8 μ L sample contain genomic DNA/ 8 μ L of positive control
 - d. 10 μ L rtPCR reaction mix

Recommended concentration: Between 0.3 to 5.0 ng/ μ L

Calculations: Mt copy number =


$\Delta Ct1 = - Ct (\text{mitochondria -control}) - Ct (\text{nucleus -control})$

$\Delta Ct2 = Ct \text{ (mitochondria - experimental)} - Ct \text{ (nucleus - experimental)}$

$\Delta\Delta Ct = \text{Sample } \Delta Ct - \text{Average } \Delta Ct2 \text{ control.}$

mtDNA level change = $2^{-\Delta\Delta Ct}$


Total DNA isolated from rat NRK-52E cells



Suggested assay plate layout: n = nucleus; mt = mitochondria; BLK = blank

	1	2	3	4	5	6	7	8	9	10	11	12
A	nBLK	nS3	nS7	nS11	nS15	nS19	mtBLK	mtS3	mtS7	mtS11	mtS15	mtS19
B	nBLK	nS3	nS7	nS11	nS15	nS19	mtBLK	mtS3	mtS7	mtS11	mtS15	mtS19
C	nPC	nS4	nS8	nS12	nS16	nS20	mtPC	mtS4	mtS8	mtS12	mtS16	mtS20
D	nPC	nS4	nS8	nS12	nS16	nS20	mtPC	mtS4	mtS8	mtS12	mtS16	mtS20
E	nS1	nS5	nS9	nS13	nS17	nS21	mtS1	mtS5	mtS9	mtS13	mtS17	mtS21
F	nS1	nS5	nS9	nS13	nS17	nS21	mtS1	mtS5	mtS9	mtS13	mtS17	mtS21
G	nS2	nS6	nS10	nS14	nS18	nS22	mtS2	mtS6	mtS10	mtS14	mtS18	mtS22
H	nS2	nS6	nS10	nS14	nS18	nS22	mtS2	mtS6	mtS10	mtS14	mtS18	mtS22

Plot of C_T versus DNA concentration



References

- 1- Wein et al., *Oncology Letters*; 6: 1098-1102, 2013
2727 Second Ave. Suite 4113 – Detroit, MI 48201 – tel. 313.961.1606 – fax. 313.963.7130 – info@detroitrandd.com – www.detroitRandD.com

- 2- Santos et al., Free Radical Biology & Medicine; 51: 1849–1860, 2011
- 3- Santos et al., Invest Ophthalmol Vis Sci 2011 Nov 11; 52 (12): 8791-8798
- 4- Edwards et al., Diabetologia; 53: 160–169, 2010.

Ref 3- Santos et al Invest Ophthalmol Vis Sci. 2011 Nov 11;52 (12):8791-8798.

