## Rat Mitochondrial DNA Copy Number

### Assay Kit

(44 reactions)

Catalog Number: MCN 2

Store at -20°C.

FOR RESEARCH USE ONLY

**Introduction:** This DNA analysis kit is for the determination of rat mitochondrial DNA copy number, *in vivo* and *in vitro*, by the comparison of mitochondrial (mt) and nuclear (n) DNA measured by real-time PCR.

#### Kit Contents:

- 96 well PCR plate
- rtPCR reaction mix.
- Validated primers to quantify mitochondrial DNA (mtDNA).
- Validated primers to quantify nuclear DNA (nDNA).
- Positive control [1.825 ng/  $\mu$ l] (isolated total DNA from rat NRK-52E cells).

#### Not Included in Kit:

- DNA isolation Kit
- Nuclease-free water
- PCR Tubes and Caps

#### Thermal cycler program:

- Preprogram PCR machine for this profile:
  - a. 95°C, 10 min**(40 Cycles)**
  - b. 95°C, 15 sec
  - c. 60°C, 60 sec

**Real time PCR procedure:** The following procedure is for each 20  $\mu$ L reaction. Increase all amounts proportionally according to the total number of tubes.

- Per PCR tube (20  $\mu$ L Rx), mix the following:
  - a. 1  $\mu$ L forward primer
  - b. 1  $\mu$ L reverse primer
  - c. 8  $\mu$ L sample contain genomic DNA/ 8  $\mu$ L of positive control
  - d. 10  $\mu$ L rtPCR reaction mix

**Recommended concentration:** Between 0.3 to 5.0 ng/ $\mu$ L

**Calculations: Mt copy number =**


$\Delta Ct1 = Ct (\text{mitochondria -control}) - Ct (\text{nucleus -control})$

$\Delta Ct2 = Ct(\text{mitochondria - experimental}) - Ct(\text{nucleus - experimental})$

$\Delta\Delta Ct = \text{Sample } \Delta Ct - \text{Average } \Delta Ct2 \text{ control.}$

mtDNA level change =  $2^{-\Delta\Delta Ct}$


Total DNA isolated from rat NRK-52E cells

Suggested assay plate layout: n = nucleus; mt = mitochondria; BLK = blank

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Plot of  $C_T$  versus DNA concentration

## References

1. Mishra JS, More AS, Kumar S. Elevated androgen levels induce hyperinsulinemia through increase in Ins1 transcription in pancreatic beta cells in female rats. *Biol Reprod.* 2018 Apr 1;98(4):520-531. doi: 10.1093/biolre/i0y017. PMID: 29365042.  
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6279097/pdf/i0y017.pdf>
2. Dos Santos JM, de Oliveira DS, Moreli ML, Benite-Ribeiro SA. The role of mitochondrial DNA damage at skeletal muscle oxidative stress on the development of type 2 diabetes. *Mol Cell Biochem.* 2018 Dec;449(1-2):251-255. doi: 10.1007/s11010-018-3361-5. PMID: 29679277.  
<https://pubmed.ncbi.nlm.nih.gov/29679277/>
3. Mishra JS, Blesson CS, Kumar S. Testosterone Decreases Placental Mitochondrial Content and Cellular Bioenergetics. *Biology (Basel).* 2020 Jul 20;9(7):176. doi: 10.3390/biology9070176. PMID: 32698476.  
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC7407169/pdf/biology-09-00176.pdf>